NG-Test hCG WB

Brzi test za kvalitativno određivanje humanog horionskog gonadotropina (hCG) u punoj krvi
Samo za profesionalnu in vitro dijagnostiku

Ref: NGB-hCG-W23-017

Primena

NG-Test hCG WB je brzi vizuelni imunoesej za kvalitativno paušalno određivanje humanog horionskog gonadotropina u uzorku pune krvi koji služi kao pomoć u ranom otkrivanju trudnoće. Osetljivost NG-Test hCG WB je 10mIU/mL

Kratak pregled i objašnjenje

Humani horionski gonadotropin (hCG) je glikoprotein koji tkivo posteljice luči tokom trudnoće. Tokom normalne trudnoće, hCG se može detektovati u uzorku pune krvi već 7 do 10 dana nakon začeća. Nivo hCG nastavlja brzo da raste, često prelazeći 100 mIU/mL i predstavlja odličan marker za potvrđivanje trudnoće. Test koristi monoklonska antitela za hCG.

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Literatura:


- Batzer FR. Fertil Steril. Hormonal evaluation of early pregnancy. 1980 Jul; 34(1): 1-13.
- Catt KJ, Dufau ML, Vaitukaitis JL. Appearance of hCG in pregnancy plasma following the initiation of implantation of the blastocyst. J Clin Endocrinol Metab. 1975 Mar; 40(3): 537-40.
- Braunstein GD, Rasor J, Danzer H, Adler D, Wade ME. Serum human chorionic gonadotropin levels throughout normal pregnancy. Am J Obstet Gynecol. 1976 Nov 15; 126(6): 678-81.
- Lenton EA, Neal LM, Sulaiman R. Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. Fertil Steril. 1982 Jun; 37(6): 773-8.
- Engvall E. Enzyme immunoassay ELISA and EMIT. Methods Enzymol. 1980; 70(A): 419-39.
- Uotila M, Ruoslahti E, Engvall E. Two-site sandwich enzyme immunoassay with monoclonal antibodies to human alpha-fetoprotein. J Immunol Methods. 1981; 42(1): 11-5.
- Steier JA, Bergsjø P, Myking OL. Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion, and removed ectopic pregnancy. Obstet Gynecol. 1984 Sep; 64(3): 391-4.
- Dawood MY, Saxena BB, Landesman R. Human chorionic gonadotropin and its subunits in hydatidiform mole and choriocarcinoma. Obstet Gynecol. 1977 Aug; 50(2): 172-81.
- Braunstein GD, Vaitukaitis JL, Carbone PP, Ross GT. Ectopic production of human chorionic gonadotropin by neoplasms. Ann Intern Med. 1973 Jan; 78(1): 39-45.

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